

Technical Datasheet

VWR Latex Examination Gloves Disposable, Powder Free, Non-Sterile

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force

Cat. No.	Size	Packing	Case
112-5397	XS	100/pack, 10packs/case	1000
112-5398	S	100/pack, 10packs/case	1000
112-5399	M	100/pack, 10packs/case	1000
112-5400	L	100/pack, 10packs/case	1000
112-5401	XL	90/pack, 10packs/case	900

MATERIAL: Natural Rubber Latex

CHARACTERISTICS:

- Examination and disposable protective gloves
- Ambidextrous, textured and beaded cuff
- Single used only
- Polymer coated inner surface, powder free and non-sterile
- Color: Natural white

DIMENSION:

Glove size	Palm Width per EN455-2	Glove Length per EN 455-2	Minimum thickness
XS	≤ 80 mm	Median ≥ 240 mm	at finger = 0.12 mm at palm = 0.10 mm at cuff = 0.06 mm
S	80 +/- 10 mm		
M	95 +/- 10 mm		
L	110 +/- 10 mm		
XL	≥ 110 mm		

PHYSICAL PROPERTIES:

Force at Break	Median ≥ 6 N (per EN455-2)
Tensile Strength	Min. 18 MPa (per ASTM D3578)
Elongation at Break	Min. 650% (per ASTM D3578)

BARRIER PERFORMANCE: AQL ≤ 1.5

Technical Datasheet

POWDER CONTENT: Max. 2.0 mg/glove (per EN455-3)

PRODUCT AND QUALITY COMPLIANCE:

- EN 455, EN 420 and EN ISO 374
- Food Contact Materials Regulations (EC) 1935/2004
- Manufactured under ISO 9001, ISO 13485 and ISO 14001 quality system.

RECOMMENDED STORAGE:

- Store in original packaging in a dry and dark place
- Recommended storage temperature: 10 °C – 30 °C
- Protect from ozone. Shield open boxes from direct sun and fluorescent lighting.

PRODUCT SHELF LIFE: 3 years

COUNTRY OF ORIGIN: Malaysia

CAUTION:

- Check suitability for the intended use and damage before use, do not use damaged gloves.
- The product contains natural rubber latex, which can cause allergic reaction including anaphylactic reactions.
- The product contains dithiocarbamate and zinc mercaptobenzothiazol – not to be used in the case of hypersensitivity to these substances.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **March 2021**