



Certificate of Conformity

VWR BEAKER WITHOUT HANDLE, PMP

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Capacity (mL)	Graduations (mL)	Sub divisions (ml)	Weight (g)	Tolerance (mL)	Maximum Top Outer Diam (mm)	Height (mm)
213-0510	25	5					
213-0482	50	10	5	11.6 ± 2	± 2.5	53.8 ± 1	58.4 ± 1
213-0483	100	20	10	18.7 ± 2	± 5	70.2 ± 1	72.5 ± 1
213-0484	250	50	25	41.6 ± 2	± 12.5	90.4 ± 1	92.1 ± 1
213-0485	500	100	50	75.2 ± 4	± 25	117.2 ± 1	132 ± 1
213-0486	1000	100	50	127 ± 5	± 50	140.5 ± 2	149.2 ± 1
213-0487	2000	200	100	200.2 ± 10	± 100	171.5 ± 2	190.2 ± 2
213-0488	5000	1000	500	341 ± 20	± 250	220.1 ± 2	260.5 ± 2

Quality System Compliance: Products are manufactured, inspected and found to be in compliance with product and quality specification requirements as documented in the ISO 9001:2008 & ISO13485:2003 Quality Management System. Inspection records are reviewed and signed for product release.

QC Testing: Representative production samples are collected and inspected in accordance with current applicable standard protocols.

Test Name	Test Reference	Accuracy	Result
Autoclaving	TEST/001	Special Insp. Level-S-1, AQL -1.5	Autoclacable
Dimensional	QSF-QC-02	Special Insp. Level-1, AQL -1.5	Within Tolerance

Material: Polymethylpentene (PMP), 100 % Virgin, FDA 21 CFR Compliance, Food grade.



Certificate of Conformity

BSE/TSE Statement: The above mentioned product is manufactured without the use of Raw materials of animal origin. BSE/TSE issue does not concern this product to the best of our knowledge.

Non-Pyrogenic: Not applicable.

Dnase Rnase Free: Not applicable.

Cytotoxicity: Not applicable.

Latex free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Storage Conditions: Can be stored in room temperature in normal warehouse conditions.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **july 2017**