

Material No.: 16004-398

Material Description: VWR Microscope Slide, Unimark, 75x25x1mm, Lilac, 72/box, 20 boxes/cs

Lot No.: 47487

Compliance and Procedure Result

This manufacturing lot has been sampled and tested in accordance with Standard Operating Procedures and has been released for the following characteristics:

Specifications:

Width (mm): 25.2 +/- 0.2mm

Length (mm): 72.2 +/- 0.2mm

Thickness (mm): 0.95 +/- 1.05mm

Material:

Soda lime glass

Pre-cleaned, ready for use, autoclavable, ground edges

Lilac marking area of approx. 20mm on one side

Adhesive

Compliance:

Standard: DIN ISO 8037-1

Country of Origin: Germany

Signed on behalf of VWR By: James Osborn

Quality Employee

5/22/2023