

S.C.A.T. Europe Overview

Range of Canisters



**THE System for
Solvent Safety**



Product Group 1: Electrically Conducting Canisters

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
A	1	S 50	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	108 421
B	2	S 60 / 61	10 litres	PE-HD-el	black	Yes	UN-Y approval	107 953
C	3	S 60 / 61	10 litres	PE-HD-el	black	No	with red viewing strip	108 216
D	4	S 60 / 61	10 litres	PE-HD-el	black	No	with green viewing strip	108 217
E	5	S 60 / 61	10 litres	PE-HD-el	black	No	with yellow viewing strip	108 215
F	6	S 60 / 61	10 litres	PE-HD-el	black	No	with blue viewing strip	108 214
G	7	S 60 / 61	10 litres	PE-HD-el	black	No	with float	108 042
H	8	S 60 / 61	20 litres	PE-HD-el	black	Yes	UN-Y approval	108 027
I	9	S 60 / 61	20 litres	PE-HD-el	black	No	with float	108 043
J	10	S 60 / 61	30 litres	PE-HD-el	black	Yes	UN-Y approval	108 192
K	11	S 60 / 61	30 litres	PE-HD-el	black	No	with float	109 193
L	12	S 70 / 71	60 litres	PE-HD-el	black	Yes	UN-Y approval	107 940
M	13	S 70 / 71	60 litres	PE-HD-el	black	No	with float	107 740
N	14	S 90	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	108 420
O	15	S 60 / 61	10 litres	PE-HD-el	black	No	with holder for sensor 108 178	108 042-S1
P	16	S 50	5 litres	PE-HD-el	black	No	designed as space-saving canister	108 317

S.C.A.T. Europe Overview

Range of Canisters



**THE System for
Solvent Safety**



Product Group 2: Non-Electrically Conducting Canisters

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
A	17	GL 45	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 952
B	18	GL 45	2.5 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 950
C	19	GL 45	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 951
D	20	S 51	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 958
E	21	S 55	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 955
F	22	S 55	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 957
G	23	S 55	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	108 175
H	24	S 55	5 litres	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	107 933
I	25	S 60 / 61	12 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 731
J	26	S 60 / 61	20 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 956
K	27	S 60 / 61	20 litres	PE-HD	natural (white)	No	UN-Y approval	108 056
L	28	S 60 / 61	30 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 959
M	29	S 60 / 61	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	108 227
N	30	S 60 / 61 & GL 45	30 litres	PE-HD	natural (white)	No	2 openings	108 115
O	31	S 65	10 litres	F-PE-HD	natural (white)	Yes	UN-Y approval	108 189
P	32	S 70 / 71	60 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 710
Q	33	S 90	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	108 020
R	34	S 95	20 litres	PE-HD	natural (white)	No		107 733
S	35	S 95	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 707


#200web / 11012021

S.C.A.T. Europe Overview

Range of Canisters



**THE System for
Solvent Safety**



107 335



107 334



107 335



107 336

Product Group 3: Space-Saving Canisters

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
A	36	S 50	5 litres	PP	natural (white)	No		107 998
B	37	S 50	5 litres	PP	natural (white)	No	with float	108 945

Product Group 4: Politainers

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
C	38	GL 38	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	107 331
D	39	GL 38	5 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 334	Box with UN-Y approval	107 330
E	40	S 60 / 61	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	107 332
F	41	S 60 / 61	20 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 336	Box with UN-Y approval	107 333

S.C.A.T. Europe Overview

Range of Canisters



**THE System for
Solvent Safety**



Product Group 5: Round Canisters

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
A	42	B 53	2 litres	PE-HD	natural (white)	No		107 711
B	43	B 83	10 litres	PP	natural (white)	No	with carrying handles	107 706
C	44	B 83	4 litres	PE-HD	natural (white)	No		107 712
D	45	B 83	50 litres	PP	natural (white)	No	with carrying handles	107 730
E	46	S 65	10 litres	PP	natural (white)	No	with carrying handles	107 706
F	47	S 65	10 litres	PE-HD	natural (white)	No	with handle and measuring scale	107 704
G	48	S 65	25 litres	PE-HD	natural (white)	No	with measuring scale	107 720
H	49	S 65	5 litres	PE-HD	natural (white)	No	with handle and measuring scale	107 722
-	50	S 65	60 litres	PE-HD	natural (white)	No	with measuring scale	107 721

Product Group 6: Flat Canisters

Image	#	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Part Number
J	51	B 63	8 litres	PP	natural (white)	No	with drain valve	107 709
K	52	S 70 / 71	20 litres	PE-HD	natural (white)	Yes	UN-Y approval	107 713