

Date of Certificate: 9/10/2025

Catalog Number: 36902-0000

Description: STERILEWARE,SCOOP,PS,2OZ

Lot Number: J397007

Date of Manufacture: 8/12/2025

Country of Origin: US

Material: High Impact Polystyrene (HIPS)

We confirm that the material manufactured with HIPS conforms to:
21 CFR 177.1640
21 CFR 178.3297

Gamma Irradiation: Irradiation exposure dose (avg.) = 23.5 min - 50.0 max kGy
ANSI/AAMI/ISO 11137 (latest revision), Method 1
No testing is done for Pyrogens / Endotoxins.

Product Shelf Life: N/A - Product should remain in its irradiated state as long as outer packaging remains completely intact.

BSE/TSE Statement: Based on the information provided to us, the HIPS resin used in the manufacture of the above product is BSE/TSE free.
During manufacture, this product does not come into contact with any products or processes that may pose a BSE/TSE risk.

Other:

Bel-Art Products hereby certifies that the above listed product has been manufactured and inspected in accordance with documented procedures and specifications, and that it conforms to all applicable requirements.

Centava Dozier
Centava Dozier
Quality Assurance Manager

R55480102

 08/29/2025 19:16:16 GMT
 Page 1 of 1

Customer Name:	Bel-Art Products-NJ	Processing Facility:	Salem	Work Order #	3820311
P.O.#	MD77165			Sales Order #	3878787
	23.5-50.0 kGy	BelArt, Gamma Treatment		Irradiation Date/Time:	08/28/2025 21:51:00 GMT
				Irradiation Cell:	CELL A

SO Line #	Qty	UOM	Customer Item Number	Customer Item Description	Customer Lot Number	Customer Load Number
101.000	60	CA	369020000(20)	STERILEWARE,SCOOP,PS,2OZ	J397007	NA
102.000	200	CA	369022002(20)	STERILEWARE,SCOOP,RED,PS,2 OZ	J395688	NA
103.000	180	CA	369300000(30)	STERILEWARE,SAMPLING SPATULA,P	J397895	NA
104.000	100	CA	369040000(20)	STERILEWARE,SCOOP,PS,4OZ	J397447	NA
105.000	60	CA	369042004(20)	STERILEWARE,SCOOP,RED,PS,4 OZ	J397009	NA
106.000	60	CA	369022010(60)	STERILEWARE,SCOOP,RED,PS,2OZ,1	J397744	NA
107.000	24	CA	136800040(24)	Flowmi Tip Strainers, 40 uM	J394329	NA
108.000	44	CA	378480001	VELVETEEN SQUARES,REPLIC,OFF W	J398338	NA
	728	CA	Total			

Quality Test Summary

Op#	Quality Test Description	Minimum Spec	Maximum Spec	Result	Pass/Fail	-----Signed By----- User	Date /Time
450.00	Minimum Dose	23.5 kGy	50.0 kGy	24.4 KGY	Pass	FWOOD	08/29/2025 11:07:11 GMT
		Reason Code Test				FRANK WOOD	
450.00	Maximum Dose	23.5 kGy	50.0 kGy	40.9 KGY	Pass	FWOOD	08/29/2025 11:07:36 GMT
		Reason Code Test				FRANK WOOD	

Sterigenics certifies that the materials listed above (as described by the Manufacturer) received the indicated doses within the precision and accuracy of the dosimetry system employed.

Electronically Signed By: KHALIL TURNER
 Reason: Work Order Completions

Date: 08/29/2025 19:15:53 GMT