

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Buffer, reference standard - pH 4.00±0.01 (25 °C)

Product Number: B5020

MDL: MFCD00130682

TEST	Specification
pH	3.99 - 4.01
Neat at 25 °C	
Traceable to NIST SRM 185 and 186	
Note	-----
Certified traceable to NIST Standards	
Material is expiry managed - actual expiration date will be indicated on Certificate of Analysis	
Specification: PRD.3.ZQ5.10000048280	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.