

Technical Datasheet

VWR WIPES RASPY

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Applicable part numbers : 115-2398, 115-2399

Attribute	Value*	Test method
Composition :		
wood pulp	51%	
polyester	49%	
Production method	spunlace apertured	
Colour	white	
Basis weight	75 g/m ²	
Absorption	760%	
Tensile strength, Dry :		ASTM D5034
Machine direction	160,1 N	
Cross direction	87,6 N	

Characteristics :

Apertured ; Slightly abrasive on one side, soft on the other ; Strong ; Resistant to acids, solvents and chlorine.
Recyclable. Food Contact Certified.

Applications :

Food & Beverage, Automotive, Industry, Printing, Laboratories, Healthcare, Janitorial.

Formats & Packaging :

Cat No	Type	Format size	Quantity	Primary pkg	Secondary pkg	Pallet
115-2398	Industrial roll	38x27cm	1 roll x 400 pcs	Plastic film	n/a	42 rolls
115-2399	Z-folded wipes	38x27cm	10 packs x 35 pcs	Plastic bag	Cardboard case	48 cases

*** Disclaimer:** The information presented here is the result of testing and other evaluations conducted during the development of the product. This information does not constitute a specification.

Date of issue: **March 2017**