

Technical Product Specification

Product

Name: HEMATOCRIT EDTA VITREX; K2EDTA 3.3mg/ml
Item no: 164213
Version: 2023-01

Physical Dimensions

Internal Diameter:	1,15 mm	Tolerance: ±	0,05 mm
External Diameter:	1,55 mm	Tolerance: ±	0,05 mm
Length:	75 mm	Tolerance: ±	0,5 mm
Nominal Volume:	75 µl	Tolerance: ±	

Material: Soda lime glass

Treatments

Impregnation: K2EDTA 3,30 mg/ml ±30%
Colour Code: Purple
Fire polishing: Light fire Polishing in colour code end

Packing

Vial: 100 pcs.
Box: 10x100 pcs. Weight/box: 0,32 Kg
Storage: To be stored in closed vials. Dry at room Temperature.

Labelling

Internal: Vitrex
External: Vitrex

Regulation

Intended use: For hematocrit testing
ISO: ISO 12772
CE: EU Regulation 2017/746 on in vitro medical devices.

