


# CH0992LP



| Eisco Labs NEXTGEN Bunsen Burner, Stabilibase Anti-Tip Design with Handle, for use with liquid propane |                                |                    |                            |                                                        |
|--------------------------------------------------------------------------------------------------------|--------------------------------|--------------------|----------------------------|--------------------------------------------------------|
| (A) Input outside Diameter                                                                             | (B) Base size                  | (C) Total Height   | (D) Adjustment Type        | (E) Additional Feature                                 |
| 0.315" (8 mm)                                                                                          | 4.25"x 4.25"<br>(108 x 108 mm) | 6.25"<br>(158.7mm) | Air flow adjustment sleeve | Candle stick easy carry handle, Brass Flame Stabilizer |