

Certificate of Conformity

VWR Catalogue Number	See List
Description	96-Well PCR Plate, Standard
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2022-02-22

Catalogue Number	Type	Well Colour	Well Volume [mL]	Coloured Lettering	Barcoded
211-0262	Standard profile, non-skirted	Clear	0.20	No	No
211-0269	Standard profile, raised well	Clear	0.20	No	No
211-0263	Standard profile, non skirted	Assorted	0.20	No	No
211-0264	Standard profile, non skirted	Blue	0.20	No	No
211-0265	Standard profile, non skirted	Green	0.20	No	No
211-0266	Standard profile, non skirted	Purple	0.20	No	No
211-0267	Standard profile, non skirted	Red	0.20	No	No
211-0268	Standard profile, non skirted	Yellow	0.20	No	No

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material Polypropylene.

Sterilization	For products that are marked sterile, they are treated with Electron-beam (E-beam) irradiation with a target of 25 kGy. The range of irradiation is verified to be between 24 and 26 kGy. The indicator on the packaging turns red when exposed to radiation >10 kGy.
ATP Assay	Not applicable.
DNase & RNase Free	The above products are free of any detectable DNase/RNase contamination.
DNA Free	Product is DNA free.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	The acceptance level for product is 0.05 EU/mL (LAL Gel Clot Method).
Latex Statement	The above products are not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.

Storage Conditions Store in room temperature.



Shelf Life 5 years after manufacture.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature



Annex

Technical Drawing

211-0269



SECTION A-A



DETAIL B

211-0262, 211-0263, 211-0264, 211-0265, 211-0266, 211-0267 and 211-0268

