

Technical Data Sheet

J.T.Baker Plate SRW Magnetic Plate SRW Magnetic - Alternatives for KingFisher™ instruments



43001-0514 (barcode version)



43001-0515 (barcode version)



43001-0516 (barcode version)

Features

- Standard SBS footprint
- Cylindrical wells enable good sample mixing
- Slightly raised rims reduce the risk of cross-contamination
- Alphanumeric code assists registration and identification of individual samples

Plate SRW Magnetic – Alternatives for KingFisher™ instruments

VWR Catalog Number	Barcode	Volume (mL)	Packed	Pk
43001-0504	-	2	10 x 6	60
43001-0505	+	NA	30 x 2	60
43001-0506	-	0.2	4 x 15	60
43001-0514	+	2	10 x 6	60
43001-0515	-	NA	30 x 2	60
43001-0516	+	0.2	4 x 15	60

ALL THIRD-PARTY TRADEMARKS REFERENCED HEREIN ARE THE EXCLUSIVE PROPERTY OF THEIR RESPECTIVE OWNERS



VWR International LLC, Radnor Corporate Center, Building One, Suite 200, 100 Matsonford Road, Radnor PA 19087, USA
VWR International Europe BV, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

Material

Plate: Polypropylene

Technical drawing and measurements

VWR Catalog Number	Barcode	Plate dimensions			Well Bottom
		Length (mm)	Width (mm)	Height (mm)	
43001-0504	-	127.76 ± 0.25	85.48 ± 0.25	44.0 ± 0.4	V
43001-0505	+	127.76 ± 0.25	85.48 ± 0.25	44.0 ± 0.4	NA
43001-0506	-	127.76 ± 0.25	85.48 ± 0.25	14.7 ± 0.4	V
43001-0514	+	127.76 ± 0.25	85.48 ± 0.25	44.0 ± 0.4	V
43001-0515	-	127.76 ± 0.25	85.48 ± 0.25	45.0 ± 0.4	NA
43001-0516	+	127.76 ± 0.25	85.48 ± 0.25	14.7 ± 0.4	V




Plate SRW 0.2 mL with V-bottom





Plate SRW 2 mL with V-bottom



Tip comb



VWR International LLC, Radnor Corporate Center, Building One, Suite 200, 100 Matsonford Road, Radnor PA 19087, USA
VWR International Europe BV, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium