



Certificate of Conformity

VWR HYGIENE MASKS, 2-PLY NON WOVEN PP

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog No: **113-0713**

Lot No: **060493**

Colour: **White**

Material: **SPP 25g/m² - White Inner PLY**
SPP 20g/m² - White Outer PLY

QC Testing:

SPP 25g/m²

Break Tensile properties(N): **MD=45.34, CD=32.89 (ASTM D5034)**

Elongation rate (%): **MD=128.38, CD=157.58 (ASTM D5034)**

SPP 20g/m²

Break Tensile properties(N): **MD=37.65, CD=21.80 (GB/T3923.1-1997)**

Elongation rate (%): **MD=118.80, CD=137.54 (GB/T3923.1-1997)**

Other Characteristics: **Latex Free**

Manufacturing Date: **24th Aug 2013**

Expiry Date: **24th Aug 2015**

This document has been produced electronically and is valid without a signature

Date of Issue: **February 2014**