

Material	20910.328
Material description	Formaldehyde 36 % (39 % w/v)
Grade	GPR RECTAPUR stabilized with about 10 % of methanol
CAS Number	50-00-0
Molecular formula	HCHO
Molecular mass	30.03

Characteristics	Specifications
Assay	35 - 38 %
Density (20/4)	1.083 - 1.093
Heavy metals (as Pb)	≤ 20 ppm
Ignition residue (SO ₄)	≤ 0.02 %
Methanol	8 - 12 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.