

Detergents &
Cleaners

Neutral

Water Bath Conditioner

Liquid

Ready-to-use

Enzymatic

Detergents

Cleaners

Ultrasonics

Liquid

Acidic
Enzymatic

Acidic
Powder
Liquid

Ultrasonics

Alkaline

Enzymatic Ready-to-use
Lens Cleaner

Detergents

Concentrated Manual

Cleaners

Concentrated

Ready-to-use

Water Bath Conditioner

Lens Cleaner

Concentrated



460 Glennie Circle
King of Prussia, PA 19406
(800) 332-6647
www.deconlabs.com



Customer Service:
800-932-5000

MANUAL CLEANERS

Contrad® 70

The Ultimate Labware Cleaner

Contrad® 70 is a liquid detergent concentrate that gets labware clean by soaking, scrubbing and ultrasonics. It rinses away totally, leaves zero residues. Phosphate-free, disposable straight to drain, Contrad® 70 is perfect for all critical cleaning, from labware to process equipment, ultrasonic cleaning and radioactive decontamination. The original and still the best... use it and see!



Dri-Clean®

Powdered Detergent

Dri-Clean® is a powdered detergent for soak, manual, and ultrasonic cleaning applications. User friendly bottle provides easy dispensing and eliminates caking. Completely soluble, free-rinsing, and residue-free. Versatile formula removes a wide variety of tenacious contaminants.



Neutrad®

Neutral pH Liquid Cleaner

Neutrad® is a liquid detergent concentrate for hand scrubbing your labware which offers high-performance cleaning with the safety of a near-neutral pH. Phosphate-free, totally rinsable, disposable straight to drain, Neutrad® is the perfect replacement for old-fashioned caustic soap powders.



Citrad®

Acidic Liquid Detergent

Citrad® is a citric acid-based detergent concentrate ideal for removing salts, scale, and metal-oxide. Also used for metal brightening and pharmaceutical-process equipment. Completely phosphate-free, totally free-rinsing and biodegradable.



Contrex® AL

Alkaline Liquid Detergent

Contrex® AL is a concentrated alkaline liquid detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex® AL is free-rinsing, phosphate-free and disposable straight to drain.



Contrex® AP

Alkaline Powdered Detergent

Contrex® AP is a concentrated powdered labware detergent formulated for use in manual, soak or ultrasonic cleaning of labware, scientific apparatus and process equipment. Contrex® AP is completely soluble, free-rinsing and safe to dispose straight to drain. The user friendly bottle allows for easy dispensing and eliminates caking.



DeSCAL®

Acidic Neutralizing Rinse (Liquid Descaler)

DeSCAL® is an inorganic acid-based detergent and descaler concentrate for removing hard water mineral deposits, rust and crusty build up from labware, water baths and stills. Low-foaming formula allows it to be used for automatic washing applications.



Contrex® CA

Acidic Liquid Detergent

Contrex® CA is a concentrated acidic liquid detergent made for use in descaling, mineral removal, oxide removal or in applications where an acidic detergent is indicated. Formulated with organic acid, Contrex® CA can be used as a manual, soak or ultrasonic cleaner.



AUTO WASHERS

Dri-Contrad®

Powdered Autowasher Detergent

Dri-Contrad® is a powdered detergent that supercharges your automatic labware washer. Non-foaming, instantly soluble, rinsable, and chlorinated for biological control. Dri-Contrad is the only powder you need for your labware washer.



Contrad® NF

Liquid Non-Foaming Autowasher Detergent

Contrad® NF is a non-foaming super-concentrated liquid detergent formulated for automatic labware washers. Phosphate-free, and totally free-rinsing, Contrad NF maximizes the efficiency of your labware machine.



Contrex® AW

Powdered Autowasher Detergent

Contrex® AW is a concentrated low-foaming powdered cleanser formulated for use in all automatic labware washing machines. Contrex AW is completely soluble, free-rinsing and biodegradable.



Contrad® 100

Liquid CIP Cleaner

Contrad® 100 is a low-foaming detergent designed for clean-in-place systems and other high-pressure spray applications, including automatic washers. Alkaline formulation is phosphate-free, residue-free, and biodegradable.



SPECIALTY CLEANERS

ElimiNase®

RNase/DNase Remover

ElimiNase® is a non-toxic, ready-to-use solution for instantly removing RNase, DNase, unwanted DNA, and radioactive contamination from labware and equipment. Safe and effective. Eliminates need for autoclaving or using DEPC.



Optical Lens Cleaner

Optical Glass Surface Cleaner

Optical Lens Cleaner is excellent for removing oil and smudges from objectives or any optical surface. It gently cleans and will not affect a coated optic or soften the mounting that holds the objective in place. This lens cleaner is laboratory-safe and non-flammable. Regular cleaning helps extend the life of your objectives by removing surface debris that could cause scratching.



Contrex® EZ

Powdered Enzymatic Detergent

Contrex® EZ is a powdered detergent formulated for use as an enzymatic detergent for manual or soak cleaning of labware and equipment contaminated with proteins and fats. Contrex® EZ prevents corrosion, is free-rinsing and disposable straight to drain.



Jet Stream®

High Pressure Dust Remover

Jet Stream® High Pressure Dust Remover is moisture-free propelled air used to remove microscopic dust from negatives, microscope slides, lenses, balances and other laboratory equipment.



Contrex® CF

Cation / Phosphate-Free Liquid Detergent

Contrex® CF is a concentrated liquid detergent formulated for use in cleaning of equipment and electronics that are sensitive to metal cations and phosphates. Contrex® CF is free-rinsing, phosphate-free and disposable straight to drain. Can be used for soak, scrub, sonication and auto washers.



BathBrite™

Water Conditioner for Water Baths

Bath Brite™ can be used to condition the water in water baths, incubators, and other temperature controlled lab equipment. 10 to 15 drops per gallon of water will keep water clear and fresh. Using Bath Brite™ will extend the time between water replacement.



Detergents & Cleaners

additional information



= Liquid



= Powder

DeSCAL®



Contrex® CA



Citrad®



Neutrad®



Contrex® AL



Dri-Clean®



Contrex® AP



Contrad® 70



Dri-Contrad®



Contrad® NF



Contrex® AW



Contrad® 100



pH: ACIDIC

NEUTRAL

ALKALINE

SOAK, MANUAL SCRUB, ULTRASONIC CLEANING

ALKALINE

AUTO WASHERS

PRODUCT GUIDE

PRODUCT	SIZE	VWR CAT. #	CASE SIZE
Contrad 70® <i>Ultimate Labware Liquid Detergent</i>	1 Liter 5 Liters 20 Liters 210 Liters	89233-152 42000-000 89233-154 89233-156	12 x 1 L 4 x 5 L — —
Dri-Clean® <i>Powdered Detergent</i>	2 kG 50 kG	89234-802 89234-800	6 x 2 kG —
Neutrad® <i>Neutral pH Liquid Cleaner</i>	950 mL 5 Liters 20 Liters	89234-808 89234-804 89234-806	12 x 950 mL 4 x 5 L —
Citrad® <i>Acidic Liquid Detergent</i>	3.8 Liters	89234-814	4 x 3.8 Liters
Contrex® AL <i>Alkaline Liquid Detergent</i>	1 qt. 3.8 Liters	89234-824 89234-822	12 x 1 qt. 4 x 3.8 Liters
Contrex® AP <i>Alkaline Powdered Detergent</i>	4 lbs. 25 lbs.	89234-818 89234-820	9 x 4 lbs. —
DeSCAL® <i>Acidic Neutralizing Rinse (Liquid Descaler)</i>	1 Liter 3.8 Liters 115 Liters	89233-170 89233-172 89233-176	12 x 1 L 4 x 3.8 Liters —
Contrex® CA <i>Acidic Liquid Detergent</i>	3.8 Liters	89234-834	4 x 3.8 Liters
Dri-Contrad® <i>Powdered Autowasher Detergent</i>	2 kG 5 kG 12.5 kG 50 kG	89234-788 89234-792 89234-794 89234-790	6 x 2 kG 4 x 5 kG — —
Contrad® NF <i>Liquid Non-foaming Autowasher Detergent</i>	3.8 Liters 10 Liters 115 Liters 210 Liters	89233-164 89233-162 89233-166 89233-168	4 x 3.8 Liters 2 x 10 L — —
Contrex® AW <i>Powdered Autowasher Detergent</i>	4 lbs. 25 lbs. 50 lbs.	89234-830 89234-832 Decon #5650	9 x 4 lbs. — —
Contrad® 100 <i>Liquid Clean-in-place Detergent</i>	3.8 Liters 20 Liters 115 Liters 210 Liters	89233-158 89166-412 89233-160 89166-410	4 x 3.8 Liters — — —
ELIMINase® <i>RNase/DNase Remover</i>	250 mL 950 mL	10798-804 10498-216	12 x 250 mL 12 x 950 mL
Contrex® EZ <i>Powdered Enzymatic Detergent</i>	4 lbs. 25 lbs.	89234-826 89234-828	9 x 4 lbs. —
Contrex® CF <i>Cation / Phosphate-Free Liquid Detergent</i>	3.8 Liters	Decon #5801	4 x 3.8 Liters
Jet Stream® <i>High Pressure Dust Remover</i>	11 oz.	89234-812	6 x 11 oz.
Optical lens cleaner™ <i>Optical Glass Surface Cleaner</i>	6 x 2 oz.	10611-616	8 x 6/pk
Bath Brite™ <i>Water Conditioner for Water baths</i>	8 oz.	10820-596	6 x 8 oz.