

Product Specification

Meets Reagent Specifications for testing Ph. Eur. monographs

Test	Specification
Assay	97.0–100.5 %
Appearance of solution	passes test
Carbonate (as Na ₂ CO ₃)	<= 2.0 %
Chlorides (as Cl)	<= 50 ppm
Heavy Metals (as Pb)	<= 20 ppm
Identification	passes test
Iron (Fe)	<= 10 ppm
Sulfates (as SO ₄)	<= 50 ppm
Additional Specification(s):	
Potassium (K)	passes test

For Laboratory, Research or Manufacturing use.

Storage conditions: Stored in an airtight non-metallic container.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000