


Recombinant Human PCNA Protein (Functional) (A476)

Specifications:

Name:	Recombinant Human PCNA Protein (Functional)
Applications:	WB, DB, SDS-PAGE, ELISA, Functional Studies
Expression System:	Escherichia coli
Nature:	Recombinant
Protein Species:	Human
Protein Length:	Full length protein
Molecular Weight:	29 kDa
Biologically Active:	Yes
Purity:	> 98% (by SDS-PAGE).
Product Form:	Liquid
Concentration:	1.0 mg/ml
Formulation:	Supplied in 25 mM HEPES, pH 7.9, with 1 mM EDTA, 0.01% NP40, 1 mM DTT, 2 µg/ml Leupeptin, 0.1 mM PMSF, 75 mM NaCl, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Recombinant Human PCNA Protein (Functional)

Recombinant Human PCNA Protein (Functional) (A476)

Images continued:

