

Material	66368.180
Material description	Polysorbate 20
Grade	GPR RECTAPUR
CAS Number	9005-64-5
Molecular formula	C ₅₈ H ₁₁₄ O ₂₆
Molecular mass	1,227.53

Characteristics	Specifications
Appearance:	.
Yellow to brown-yellow liquid	Passes test
IR Spectrum	Passes test
Density (20/4)	1.095 - 1.105
Hydroxyl value	96 - 113
Saponification value	40 - 51

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250
Briare; France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.