



Description: Sterile hypodermic needle for single use.

Components:

Description	Raw Materials
Hub	Polypropylene + master batch
Cannula	Stainless steel
Adhesive	Epoxy resin
Lubricant	Dimethyl silicone oil



Product Range:

Product Code	Colour	Cannula Ø		Cannula Length		Bevel	
AN*1838R1	Pink	18G	1.2 mm	1 1/2 inch	38 mm	Regular	11°
AN*1838S1	Pink	18G	1.2 mm	1 1/2 inch	38 mm	Short	17°
AN*1925R1	Ivory	19G	1.1 mm	1 inch	25 mm	Regular	11°
AN*1938R1	Ivory	19G	1.1 mm	1 1/2 inch	38 mm	Regular	11°
AN*1938S1	Ivory	19G	1.1 mm	1 1/2 inch	38 mm	Short	17°
AN*1950R1	Ivory	19G	1.1 mm	2 inch	50 mm	Regular	11°
AN*2025R1	Yellow	20G	0.9 mm	1 inch	25 mm	Regular	11°
AN*2038R1	Yellow	20G	0.9 mm	1 1/2 inch	38 mm	Regular	11°
AN*2116R1	Green	21G	0.8 mm	5/8 inch	16 mm	Regular	11°
AN*2125R1	Green	21G	0.8 mm	1 inch	25 mm	Regular	11°
AN*2138R1	Green	21G	0.8 mm	1 1/2 inch	38 mm	Regular	11°
AN*2150R1	Green	21G	0.8 mm	2 inch	50 mm	Regular	11°
AN*2232R1	Black	22G	0.7 mm	1 1/4 inch	32 mm	Regular	11°
AN*2238R1	Black	22G	0.7 mm	1 1/2 inch	38 mm	Regular	11°
AN*2316R1	Blue	23G	0.6 mm	5/8 inch	16 mm	Regular	11°
AN*2325R1	Blue	23G	0.6 mm	1 inch	25 mm	Regular	11°
AN*2332R1	Blue	23G	0.6 mm	1 1/4 inch	32 mm	Regular	11°
AN*2425R1	Purple	24G	0.55 mm	1 inch	25 mm	Regular	11°
AN*2516R1	Orange	25G	0.5 mm	5/8 inch	16 mm	Regular	11°
AN*2525R1	Orange	25G	0.5 mm	1 inch	25 mm	Regular	11°
AN*2613R1	Brown	26G	0.45 mm	1/2 inch	13 mm	Regular	11°
AN*2623R1	Brown	26G	0.45 mm	7/8 inch	23 mm	Regular	11°
AN*2716R1	Grey	27G	0.4 mm	5/8 inch	16 mm	Regular	11°
AN*2719R1	Grey	27G	0.4 mm	3/4 inch	19 mm	Regular	11°
AN*3013R1	Yellow	30G	0.3 mm	1/2 inch	13 mm	Regular	11°

Packaging:

Single unit:

Blister pack (film/paper)

Unit box:

100 pieces

Dimensions = 105 x 84 x 85 mm

For item numbers AN*1950R1 and AN*2150R1:

Dimensions = 105 x 84 x 110 mm

Shipping carton:

10.000 pieces

Dimensions = 450 x 440 x 450 mm

For item numbers AN*1950R1 and AN*2150R1:

5.000 pieces

Dimensions = 545 x 240 x 450 mm

Sterilization:

With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Storage:

Shelf life of 5 years

* The batch number and expiry date are printed on the packaging

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use

For single use only

Do not use if unit package is damaged.

Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Product Service GmbH, TÜV SUD Group 0123

Manufacturer:

Zhejiang Kindly Medical Devices Co. Ltd. (China)

EC representative:

Shanghai International Corp. GmbH (Europe) - Hamburg, Germany

Distributor:

TERUMO EUROPE N.V.

Interleuvenlaan 40

B - 3001 Leuven



Terumo Corporation

44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.

Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

European Offices

Terumo Europe N.V.

Benelux Sales Division
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.

Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex, France
Tel.: +33 130 96 13 00
Fax: +33 130 43 60 85

Terumo Deutschland GmbH

Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH

Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Europe España SL

Avda. Juan Carlos I, Nº 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares (Madrid),
Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Italia SRL

Via Simone Martini, 143 /145
00142 Roma RM,
Italy
Tel.: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo Sweden AB

Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
426 71 Västra Frölunda,
Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

Terumo Denmark

Filial of Terumo Sweden AB
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

Terumo UK Ltd

Tamesis
The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 220
Fax: +44 1784 476 201

Terumo Europe N.V.

Moscow Representative Office

BC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: + 7 495 988 4740
Fax: + 7 495 988 4739

Terumo Corporation Dubai Branch

(Middle East and Africa)

Al Masraf Tower, 22nd Floor
P. O. Box 20291 Dubai
U.A.E.
Tel.: + 971 4 221 2220
Fax: + 971 4 221 3330

**Addresses subject to change,
please consult our Website.**

© Registered Trademark
Published by Terumo Europe N.V.
HO 131GB-0912TE-I(09.12)E