



Product Specification Sheet

VWR CRITICAL SWAB, WITH COTTON HEADS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Description	Handle length (mm)	Handle width (mm)	Head length (mm)	Head width (mm)	Total length (mm)	Pk
149-0331	Absorbent cotton tip and wooden shaft, autoclavable	152	2.5	15.9	4.8	168	1.000
149-0332	Large absorbent cotton tip and wooden shaft, autoclavable	152	2.5	17.4	5.9	170	500
149-0340	Double cone-shaped, highly absorbent cotton heads, packaged in an anti-static leaf pack, white	82.0	2.5	11	4.7	93	1.250
149-0341	Double, highly absorbent cotton heads, packaged in an anti-static leaf pack, white	78	2.5	15	4.4	93	2.500
149-0342	Double slim highly absorbent cotton heads, packaged in an anti-static leaf pack, white	75.5	1.5	14	3.2	89.5	2.500





Product Specification Sheet

VWR CRITICAL SWAB, WITH COTTON HEADS

Applications:

- Cleaning disk drives, optics, PC board assembly and other critical industrial components
- Cleaning sensitive components
- Veterinary
- Weapons cleaning

Other Characteristics:

- Highly absorbent
- Low-linting formulations reduce contamination
- Economical

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.

Date of issue: **April 2013**



VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Version: 2.0