

Corning® Stripwell™ Strips and Microplates to Support Serological Testing Workflows


Detection of antibodies to infectious agents such as the SARS-CoV-2 virus that causes COVID-19 is routinely performed using the enzyme-linked immunosorbent assay (ELISA). Numerous serological test kits for COVID-19 have been developed and received emergency use authorization (EUA) from regulatory agencies in the past few months. Many of these test kits utilize ELISA assays to identify the presence of antibodies to the SARS-CoV-2 virus in blood samples.

Corning is an industry leader in the breadth of high quality, high-performance ELISA 1 x 8 strips and microplates available to support laboratory testing workflows. These ELISA products can be used as components in a wide range of testing workflows. Surface chemistries are also available for ELISA workflows requiring covalent immobilization to sulfhydryl groups or carbohydrate moieties.

ELISA strips and microplates are available in clear, black, or white polystyrene with a variety of different surface chemistries to suit multiple assay formats. In addition, Corning has expanded manufacturing capabilities of all ELISA products to ensure global availability of these critical components in diagnostic workflows.

Features and Benefits

- ▶ Industry-leading manufacturing capabilities for high quality and high-performance ELISA strips and microplates
- ▶ Offered in a wide variety of surface chemistries to cover a multitude of assay formats
- ▶ High binding and medium binding surfaces
- ▶ Available manufacturing capacity to support global laboratory supplies requirements for COVID-19



Ordering Information

VWR Cat. No.	Corning Cat. No.	Description	Qty/Pk	Qty/Cs
29442-302	2592	96-well clear polystyrene microplate, 1 x 8 Corning® Stripwell™ strips, high bind, without lid, nonsterile	25	100
29442-304	2593	96-well clear polystyrene microplate, 1 x 8 Corning Stripwell strips, medium bind, without lid, nonsterile	25	100
29444-024	3913	96-well white polystyrene microplate, 1 x 8 Corning Stripwell strips, medium bind, without lid, nonsterile	25	100
29444-028	3923	96-well white polystyrene microplate, 1 x 8 Corning Stripwell strips, high bind, without lid, nonsterile	25	100
29444-026	3914	96-well black polystyrene microplate, 1 x 8 Corning Stripwell strips, medium bind, without lid, nonsterile	25	100
29444-030	3924	96-well black polystyrene microplate, 1 x 8 Corning Stripwell strips, high bind, without lid, nonsterile	25	100
46610-186	2481	96-well clear polystyrene microplate, 1 x 8 Corning Stripwell strips, high bind, low volume, without lid, nonsterile	25	100
10068-004	2580	96-well clear polystyrene microplate, 1 x 8 Corning Stripwell strips, high bind, without frame or lid, nonsterile, bulk packaged	--	800
29442-282	2572	96-well Corning Stripwell egg crate (microplate) strips holder	--	20

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.




Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by country, federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Offers valid in countries listed above, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.