

Material	22355.367
Material description	Calcium hydroxide
Grade	GPR RECTAPUR
CAS Number	1305-62-0
Molecular formula	Ca(OH) ₂
Molecular mass	74.09

Characteristics	Specifications
Assay (at filling)	≥ 96 %
Heavy metals (as Pb)	≤ 50 ppm
Insolubility in hydrochloric acid	≤ 0.1 %
Not precipitated by (COONH ₄) ₂ (as SO ₄)	≤ 2.5 %
Cl (Chloride)	≤ 100 ppm
CO ₃ (as CaCO ₃)	≤ 3.0 %
SO ₄ (Sulphate)	≤ 0.2 %
Fe (Iron)	≤ 0.05 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.