

Certificate of Conformity

VWR Catalogue Number	See List
Description	PHP Plates, Sterile
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-07-08

VWR Catalogue Number	Description	Diameter [mm]	Pk
VWRA391-0645	PHP petri plate, stack size 10	90	600
VWRA391-0662	PHP petri plate, stack size 10, double wrapped	90	600
VWRA391-0685	PHP petri plate, stack size 10	90	600
VWRA391-0686	PHP petri plate, stack size 10	90	600
VWRA391-0687	PHP petri plate, stack size 10	90	600
VWRA391-0691	PHP petri plate, stack size 10	90	600
VWRA391-0692	PHP petri plate, stack size 10	90	600
VWRA391-0684	PHP petri plate, stack size 10	55	1.120
VWRA391-0646	PHP petri plate, stack size 20	90	600
VWRA391-0663	PHP petri plate, stack size 20, double wrapped	90	600
VWRA391-0688	PHP petri plate, stack size 20	90	600
VWRA391-0689	PHP petri plate, stack size 20	90	600
VWRA391-0693	PHP petri plate, stack size 20	90	600
VWRA391-0694	PHP petri plate, stack size 20	90	600

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is manufactured in Class 100K (ISO 7) clean room, compliant with ISO 14644-1 standards.

Material	Virgin Crystal Polystyrene; US FDA approved		
Dimensional Specification	Lid Diameter	91.20 mm	
	Lid Height	7.50 mm	
	Dish Diameter	85.30 mm	
	Dish Height	12.18 mm	
Sterilization	Products are Gamma irradiated as per ISO 11137 which meets a minimum requirement of 10 ⁻⁶ SAL.		
ATP Assay	Not Applicable.		
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.		
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).		
Cytotoxicity	This product is free of cytotoxicity effects.		
Non-Pyrogenic Statement	The acceptance level for product is ≤ 0.01 EU/mL (LAL Gel Clot Method).		
Latex Statement	This product is not made from natural rubber latex.		
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.		
DEHP Statement	Not Applicable.		
RoHS	Product comply with the limits as set by RoHS Directive (EU) 2015/863 amending Annex II to Directive 2011/65/EU.		
REACH Statement	Product meet the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006. To the best of our knowledge product does not contain SVHC's above the regulatory threshold of 0.1%.		

Storage Conditions Room temperature.

Shelf Life 3 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	07-2025
2	2.0	Modified Quality system compliance and disclaimer	07-2025