



Sol - Care® Safety Multi-Sample Needle

- Safety shield aligns with a bevel up needle orientation, allowing for low-angle penetration
- Pivoting safety shield with audible click that indicates it is locked in place
- Safety mechanism is designed for one-handed activation
- Safety shield is color coded to reduce potential errors

***samples available upon request**

VWR Cat No. Product Description

76582-902 Sol-Care Safety Multi-Sample Needle 22G x 1-1/4"

76582-900 Sol-Care Safety Multi-Sample Needle 21G x 1-1/4"

VWR Cat No. Product Description

76582-906 Sol-Care Safety Multi-Sample Needle with Pre-attached Holder 22G x 1-1/4"

76582-904 Sol-Care Safety Multi-Sample Needle with Pre-attached Holder 21G x 1-1/4"