

Certificate of Conformity

VWR Catalogue Number	115-0049
Description	Multi-layer adhesive floor mats, PureStep 30 layers, 8 µm ±20% mm thickness, 610 mm Length, 762 mm Width, Grey
Lot Number	UT0304
Expiry Date (yyyy-mm)	2029-08
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-07-21

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 14000:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Low Density Polyethylene, (LDPE) Blend, with 20% biohybrid resin.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.

Non-Pyrogenic Statement	Not Tested.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store at 12-35 °C (54-95F).
Shelf Life	Minimum 5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Annex

Mat Layers: Low Density Polyethylene.

1.8 mil (45 micron) for 30 Layer, 4 mat**

1.4 mil (35 micron) for 30 Layer, 8 mat and 60 Layer**

Adhesive: 0.3 mil (7 micron) acrylic based**

Peel Adhesive: 9.0 oz. (250g)/25mm, Individual layers (+/- 50g)

24 oz. (700g)/25mm, Bottom layer. (+/- 50g)

Tensile Strength: >1.2 kg/cm

Heat Resistance: 158° F/70° C -48 hours

Surface		Bottom	
Tack Spec	Tack Actual	Tack Spec	Tack Actual
150-230	186	550-650	614

Other Characteristics:

- Anti-microbial, destroys germs and fungi
- Non-slip, colour-coded
- Reduced maintenance and cleaning costs