

# Certificate of Conformity

|                                           |                                                        |
|-------------------------------------------|--------------------------------------------------------|
| <b>VWR Catalogue Number</b>               | 732-1614                                               |
| <b>Description</b>                        | Filter tips, Low retention, Extra long, 10 µl, Sterile |
| <b>Lot Number</b>                         | 40852-008C4-008B                                       |
| <b>Date of Manufacture (yyyy-mm-dd)</b>   | 2020-02-16                                             |
| <b>Date of Sterilization (yyyy-mm-dd)</b> | 2020-02-21                                             |
| <b>Expiry Date (yyyy-mm-dd)</b>           | 2023-02                                                |
| <b>Country of Origin</b>                  | Manufactured in USA                                    |
| <b>Date of Issue (yyyy-mm-dd)</b>         | 2021-11-23                                             |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                               |                                                                                                                                                                                                                                                                                   |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | 100% Virgin Polypropylene.<br><br>Materials used meet FDA Regulations 21CFR Compliance. Resins are USP Class VI certified.                                                                                                                                                        |
| <b>Sterilization</b>          | Product labeled as sterile is gamma or beta irradiated and dose released upon ANSI/AAMI/ISO 11137 recommended practices in effect at the time of validation. Products labeled sterile meets a minimum requirement of $10^{-6}$ SAL with a calculated dose range of 15.9–27.5 kGy. |
| <b>ATP Assay</b>              | This product has been tested and is ATP free. The assay detection limit is less than $2 \times 10^{-12}$ mg/µL.                                                                                                                                                                   |
| <b>DNase &amp; RNase Free</b> | This product is free of any detectable DNase/RNase contamination, Human DNA and Protease.                                                                                                                                                                                         |

|                                |                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).          |
| <b>Cytotoxicity</b>            | Testing is conducted and have been shown to be non-toxic.                                                                                                                                                                                                                                                                                                            |
| <b>Non-Pyrogenic Statement</b> | The acceptance level for product is less than 0.05 EU/mL (LAL Assay).                                                                                                                                                                                                                                                                                                |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                                  |
| <b>BPA Statement</b>           | Bisphenol are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                    |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                   |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed. |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.                                                                                                                                                      |
| <b>Storage Conditions</b>      | Store at room temperature.                                                                                                                                                                                                                                                                                                                                           |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.