

# Certificate of Quality

|                                   |                                           |
|-----------------------------------|-------------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                                  |
| <b>Description</b>                | VWR® Microscope Slide Boxes for 50 Slides |
| <b>Country of Origin</b>          | Manufactured in China                     |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2026-03-26                                |

| VWR Catalog Number | Color | Unit       |
|--------------------|-------|------------|
| 82024-580          | Blue  | Each       |
|                    |       | Pack of 24 |
|                    |       | Case of 72 |
| 82024-582          | Red   | Each       |
|                    |       | Pack of 24 |
|                    |       | Case of 72 |
| 82024-584          | White | Each       |
|                    |       | Pack of 24 |
|                    |       | Case of 72 |

## Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared that the products are inspected and controlled through whole-production processing in accordance with current applicable product specifications and quality control SOPs. Inspection records are reviewed and signed off by qualified personnel for product release.

## QC Testing

Representative products samples are collected and inspected in accordance with AQL level II 4.0 sampling plan per MIL-STD-1916 and current applicable product specifications.

## Product Specifications

|                               |                                                                                                                                                                                                 |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Injection-molded of polypropylene homopolymer with non-toxic blue, red, white pigment (refer table for exact color) which are inert or minimally reactive with most common laboratory reagents. |
| <b>Sterilization</b>          | The products are non-sterile.                                                                                                                                                                   |
| <b>ATP Assay</b>              | Not applicable.                                                                                                                                                                                 |
| <b>DNase &amp; RNase Free</b> | Not applicable.                                                                                                                                                                                 |

|                                |                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE). |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                              |
| <b>Latex Statement</b>         | The products are latex free.                                                                                                                                                                                                                                                                                                                                 |
| <b>BPA Statement</b>           | Bisphenols are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                           |
| <b>DEHP Statement</b>          | Products comply with the requirements in the US Consumer Product Safety Improvement Act (CPSIA) of 2008 and the European Union Regulation (EC) 1907/2006 regarding the concentration limit and the usage of the phthalate (DEHP, CAS RN 117-81-7).                                                                                                           |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routine analysis is performed.                             |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routine analysis is performed.                                                                                                                                                                    |
| <b>Storage Conditions</b>      | Storage temperatures -20°C to 60°C (-4°F to 140°F)<br>Relative humidity during storage: 10 % to 90 %                                                                                                                                                                                                                                                         |
| <b>Shelf Life</b>              | Indefinite shelf life when stored in the recommended conditions.                                                                                                                                                                                                                                                                                             |

Disclaimer: VWR'S Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co- manufacturing, or processing activities such as design, packaging, labelling, testing, or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use. This document has been produced electronically and is valid without a signature.

### Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoQ using new layout | March 2026   |