

Certificate of Analysis

Leptin, Human, Recombinant, *E. coli*

Batch Number: 3154000
Material Number: 429700-1MG
Molecular Weight: 16000
CAS Number: 177404-21-6

Quality Release Date: 07-DEC-2018

Recommended Retest Date: 30-NOV-2023

Analytical Data

Test	Tolerance	Result
Formulation:		Lyophilized from 650.0 µL of a 0.2 µm filtered solution in 10 mM Tris-HCl and 150 mM NaCl, pH 8.5.
Form:		Solid
Purity by SDS-PAGE:	≥97.0 %	>97.0 %
Endotoxins:	≤1.0 EU/µg	<0.1 EU/µg

Storage and Handling: ≤ -70°C

This lot conforms to specifications established by EMD Millipore Corporation for this product.



Quality Control / Assurance Signature

07-DEC-2018

Date

Prices and availability are subject to change. ©Copyright 2018 EMD Millipore Corporation, an affiliate of Merck KGaA, Darmstadt, Germany. All rights reserved. Each product is sold with a limited warranty, which is provided with each purchase. Each product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes, nor is it intended for human use. EMD Millipore products may not be resold, modified for resale, or used to manufacture commercial products without written approval of EMD Millipore.

EMD Millipore Corporation | 28820 Single Oak Dr., Temecula, CA 92590

Technical Support NA +1-800-221-1975 | email: www.millipore.com/techservices |

www.calbiochem.com

Technical Support All Other Countries - Contact Your Local Office

FOR RESEARCH USE ONLY.

Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purpose other than research is strictly prohibited.

Calbiochem and all other trademarks, unless specifically identified as belonging to a third party, are owned by Merck KGaA, Darmstadt, Germany