



Digital Blood Pressure Monitor

Competitive Cross Reference

Model / Features	ADC e-sphyg 3	WA ProBP 2400	Omron 907XL
Manufacturer	ADC	WA	Omron
Measures systolic/diastolic/pulse	Y	Y	Y
Measures MAP	Y	Y	N
Static pressure accuracy	3 mmHg	3 mmHg	3 mmHg
Pulse accuracy	5%	5%	5%
Adjustable initial inflation pressure	Y	Y	Y
Microprocessor controlled deflation	Y	Y	Y
Irregular heartbeat indicator	Y	Y	Y
Average mode	Y	Y	Y
Constituent readings of avg mode	Y	Y	Y
Manual mode	Y	Y	Y
Dual mode display	N	N	N
Backlit display	Y	Y	N
Standard high-capacity battery	Y	N	N
Measurements/full charge	875-1020	600-700	300
User-enabled calibration check	Y	Y	Y
Last-reading memory	Y	Y	Y
Number of cuffs included	3	2	4
Requires specialized cuffs	N	Y	Y
Range of cuffs	19-50 cm	22-42 cm	17-50 cm
Overall tubing length	190 cm	122 cm	130 cm
Standard bayonet cuff connectors	Y	N	N
Battery type	NiMh	NiMh	NiMh
AC adapter with international plugs included	Y	Y	N
Main unit dimensions (LxWxH)	20 x 12.5 x 9 cm	20 x 12.5 x 9 cm	20.3 x 13.9 x 13.1 cm
Main unit weight, including battery	942 g	800 g	910 g
Instrument warranty	3 years	2 years	5 years
Cuff warranty	3 years	1 year	1 year
MSRP	\$429.99	\$435.00	\$749.99

Platforms

Number of optional platforms	3	2	3
Height-adjustable mobile stand	Y	N	Y
Mobile stand MSRP	\$209.99	\$278.00	\$105.00
Wall mount	Y	Y	Y
Wall mount compatible with mfg wall system	Y	Y	N

