

5 mL, Model 1005 LTN CTC SYR, Cemented NDL, 23 ga, point style 5 (110 °C)

Specification Sheet (Part/REF # 203086)



Product Specifications

Volume	5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Needle (LTN)
Needle	Cemented Needle
Gauge	23 gauge
Point Style	5
Needle Length	56 mm (2.20 inches)
Fluid Path	SST, Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	29 psig (2bar)
Maximum Temperature	230°F (110°C)
Minimum Temperature	68°F (20°C)
Needle Inner Diameter	0.008 inches (0.20 mm)
Needle Outer Diameter	0.025 inches (0.64 mm)
Std Needle Dead Volume	4.53 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/syringes/203086>

Spec. Version B