

Light Pack

The Responsible Tip Box



Eco-conscious



Reduce your lab's
environmental footprint



Forward – thinking

We've created the responsible pipette tip box: a sturdy cardboard design packed with environmental benefits – without losing performance.

- Compared to the single tray plastic rack, the plastic content of the box has been reduced by 74%
- The cardboard of the Light Pack is FSC certified
- The tray is made from recyclable plastic
- The window is made from acetate cellulose
- Awarded with the My Green Lab ACT Ecolabel

All this means that you have a pipette tip box that delivers measurable impact. When you're done, the box folds flat for efficient stacking and easy integration into established recycling workflows.



reddot winner 2026

>>

Ordering Information

Volume	Purity Level	Tips/Units	Order Code
200µl	▪ Free of DNase, RNase, Human DNA and Endotoxins	▪ Optifit ▪ 10 x 96	LH-G790200
200µl	▪ Free of DNase, RNase, Human DNA and Endotoxins ▪ Pre-sterilized	▪ Optifit ▪ 10 x 96	LH-G790201
200µl	▪ Free of DNase, RNase, Human DNA and Endotoxins ▪ Pre-sterilized	▪ Safetyspace ▪ 10 x 96	LH-G790201F
120µl	▪ Free of DNase, RNase, Human DNA and Endotoxins ▪ Pre-sterilized	▪ Safetyspace ▪ 10 x 96	LH-G790101F

Optilid
Hinged lid which can be opened and closed with one hand

Cardboard box with Plastic Tray
Plastic reduction by 74%

Cellulose Acetate Window
to see the number of tips in the box

Foldable
At end of use, reducing waste volume by up to 50%

Safety Corner and Seamless Bottom
Cleanliness and reduction in contamination risk

Robust and Rigid Design
For stability during tip picking

Germany
Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Straße 20
37079 Göttingen
Phone +49 551 308 0

USA
Sartorius Corporation
3874 Research Park Drive
Ann Arbor, MI 48108
Phone +1 734 769 1600