

Certificate of Conformity

VWR Catalogue Number	See List
Description	Test tubes without rim
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2021-10-07

Catalogue Number	Capacity [ml]	Length [mm]	External Diameter [mm]	Material	Type
216-1163	3	55	11.5	Polystyrene	Non graduated
216-1164	3	55	11.5	Polypropylene	Non graduated
212-0668	3	70	11	Polystyrene	Non graduated
216-1165	5	75	12	Polypropylene	Sorwall type
212-0728	5	75	12	Polystyrene	Non graduated Sorwall type
212-0729	5	75	12	Polystyrene	Non graduated Sorwall type
212-0727	5	75	13	Polystyrene	Non graduated
216-1167	5	75	13	Polypropylene	Non graduated
216-1044	10	100	16	Polystyrene	Non graduated
216-1046	10	100	16	Polypropylene	Graduated
216-2213	10	100	16	Polystyrene	Graduated

Stoppers

Catalogue Number	Top Diameter [mm]	Colour
212-0679	12	Natural
212-0572	12	White
212-0573	12	Blue
212-0574	12	Red
212-0575	12	Green
212-0576	12	Yellow
212-0666	13	Yellow
212-0665	13	Green
212-0664	13	Natural

Catalogue Number	Top Diameter [mm]	Colour
212-0667	13	Red
212-0577	13	Blue
212-0578	13	White
212-0550	16	Natural
212-0551	16	White
212-0552	16	Blue
212-0553	16	Green
212-0554	16	Red

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material	See above table. The products are considered safe to use in food, pharmaceutical and cosmetic contact applications.
Sterilization	Non Sterile.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.

Non-Pyrogenic Statement	Not applicable.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Bisphenol A (BPA) has not been intentionally incorporated in the manufacturing of the products.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store in a dry place, away from sources of heat and direct sunlight.
Shelf Life	60 months.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature