

Certificate of Conformity

VWR Catalogue Number	613-0230
Description	Pipette tips, 1 – 200 µl, Bevel Point (yellow) with 49.66 mm Length, Sterile
Lot Number	41048-111C4-111AJ
Date of Sterilization (yyyy-mm-dd)	2021-05-12
Expiry Date (yyyy-mm-dd)	2024-03
Country of Origin	Manufactured in USA
Date of Issue (yyyy-mm-dd)	2024-01-19

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	100% Virgin Polypropylene, Polyethylene or Polystyrene. Materials used meet FDA Regulations 21CFR Compliance.
Sterilization	Product labeled as sterile is gamma or beta irradiated within a calculated dose range of 17.1–26.9 kGy. Dose range is necessary for the stated sterility assurance level is continuously audited through quarterly bio-burden and sterility validation studies performed according to ANSI/AAMI/ISO 11137 by an independent laboratory. Products labeled sterile meets a minimum requirement of 10^{-6} SAL.
ATP Assay	Product sample surface are tested for present of ATP. The test result is less than 2×10^{-12} mg/µL of ATP for the product be certified ATP free.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Testing is conducted and have been shown to be non-toxic.
Non-Pyrogenic Statement	The acceptance level for product is less than 0.05 EU/mL (LAL Assay).
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Store at room temperature.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.