

Certificate of Quality

VWR EU Catalogue Number	VWRI210-000033
Description	VWR®, Cryo coder cap inserts, assorted
Lot Number	V19H24
Date of Manufacture (yyyy-mm)	2024-08
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-06-06

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 and ISO 13485:2016** standards. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level-1, AQL -1.5	Autoclavable
Dimensional	QWI-QC-23	General Insp. Level-1, AQL -1.5	Within Tolerance

Product Specifications

Material	Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed
Temperature Resistance	Not Applicable
Sterilization	Not Applicable
ATP Assay	Not Applicable
DNase & RNase Free	Not Applicable

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable
Non-Pyrogenic Statement	Not Applicable
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Not Applicable
DEHP Statement	Not Applicable
RoHS	Not Applicable
REACH Statement	Not Applicable
Storage Conditions	Store and transport at ambient Temperature. Do not expose to extreme temperature and direct Sunlight.
Shelf Life	Not Applicable

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.