

Certificate of Conformity

VWR Catalogue Number	See List
Description	Shell vials
Country of Origin	Manufactured in Germany - 548-0042A, 548-0043A, 548-0347A, 548-0352A, 548-0353A, 548-2481A, 548-0045A, 548-0046A, 548-0449A, 548-0048A, 548-0049A, 548-0044A, 548-0050A Manufactured in China - 548-0830A Manufactured in USA - 548-2411A, 548-1389A
Date of Issue (yyyy-mm-dd)	2022-12-29

Vials, Flat bottom base, Borosilicate glass 5.1, No cap included

Catalogue Number	Description	Capacity [mL]	Ø x H [mm]	Colour
548-0042A	Shell vials, 8 mm PE plug with insertion barrier	1.00	8 x 40	Clear
548-0043A	Shell vials, 8 mm PE plug with insertion barrier	1.00	8 x 40	Amber
548-0347A	Shell vials, 6 mm PE plug	1.00	7.8 x 35	Clear
548-0352A	Shell vials, 8 mm PE plug	1.00	8 x 40	Clear
548-0353A	Shell vials, 8 mm PE plug	1.00	8 x 40	Amber
548-2481A	Shell vials without plug	1.00	8 x 40	Clear
548-0045A	Shell vials, 12 mm PE plug	2.00	12 x 31.5	Clear
548-0046A	Shell vials, 12 mm PE plug	2.00	12 x 31.5	Amber
548-0449A	Shell vials without plug	2.00	12 x 31	Clear
548-0048A	Shell vials, 15 mm PE plug	4.00	15 x 45	Clear
548-0049A	Shell vials, 15 mm PE plug	4.00	15 x 45	Amber

Inserts, Conical bottom base, Borosilicate glass 5.1, No cap included

Catalogue Number	Description	Capacity [mL]	Ø x H [mm]	Colour
548-0044A	Micro-insert, 13 mm top	0.10	5 x 34	Clear
548-0050A	Micro-insert, 13 mm top	0.30	6 x 43.5	Clear

Plugs

Catalogue Number	Description	Cap Colour	Cap Size [mm]	Cap Material
548-2411A	15 mm PE-plug	Transparent	15	PE
548-0830A	8 mm PE-plug with insertion barrier	Transparent	8	PE
548-1389A	8 mm PE-plug	Transparent	8	PE

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Refer table.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.

Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Not applicable.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	12-2022