

Certificate of Conformity

Catalogue Number	43001-05041
Description	PLATE 96 SRW 2 ML PP CS60
LOT no.	1154304/1-1
Country of Origin	Manufactured in Germany
Date of Issue (dd/mm/yyyy)	09/01/2026
Production date	2025/12 – 2026/01
Expiry date^{10/}	2030/12 – 2031/01

Quality System Compliance

Products are manufactured in [ISO 13485:2016](#) and [ISO 9001:2015](#) certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polypropylene
Manufacturing Environment	Production lines are subject to hygiene controls. Staff wear protective clothing. Products are regularly tested for purity by an external laboratory.
Sterilization	Not applicable.
ATP Assay	Not applicable.
PCR inhibitors Free	Product yield equivalent to control < 10 amplifiable targets
Human DNA Free	Tested* with detection level < 1 pg
DNase & RNase Free	Tested* with detection level < 1*10 ⁻⁶ KunitzU (DNase) and < 1*10 ⁻⁹ KunitzU (RNase)
BSE/TSE	Tallow derived additives may be used in the manufacture of the raw materials. These tallow derived additives are in compliance with the guidance EMA/410/01 Rev. 3 – July 1, 2011 "Minimizing the risk of transmitting animal spongiform encephalopathy agents via human

and veterinary medicinal products", adopted by the European Commission and published in the Official Journal of the European Union March 5, 2011 (2011/C 73/01).

To the best of our knowledge, these treatments ensure inactivation of any agents present in the tallow derivatives to be used in plastics materials. We therefore can state that this grade is considered safe to use with respect to BSE and TSE transmissions.

Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	No intentionally use in the formulation of the product.
BPA Statement	No intentionally use in the formulation of the product.
DEHP Statement	Not applicable.
Endotoxine Free	Tested* with detection level < 0,00033 EU/Cavity
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).
REACH Statement	The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.
Material Statement	Not applicable.
Storage Conditions	Keep dry and away from sunlight.
Shelf Life	5 years.

* skip lot testing

Disclaimer: Ritter states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.