

Certificate of Quality

VWR NA Catalogue Number	89236-082
VWR EU Catalogue Number	479-2414
Description	VWR®, Cryo Boxes, Polycarbonate, 81 compartments, for 4 -5 mL tubes
Lot Number	V30L24
Date of Manufacture (yyyy-mm)	2024-12
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2025-06-16

Quality System Compliance

Products are manufactured under the [ISO 9001:2015 & ISO 13485:2016](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Test Name	Test Reference	Accuracy	Result
Autoclaving	QWI-QC-07	General Insp. Level- I, AQL- 1.5	Autoclavable
Dimensional	QWI-QC-23	General Insp. Level- I, AQL- 1.5	Within tolerance

Product Specifications

Material	Polycarbonate.
Sterilization	Non-sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	The products are not made from natural rubber latex.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Application temperature -196 °C to +121 °C. Product to be kept in a cool, dry rack at room temperature. The product should not be exposed to the sunlight or extreme temperatures.
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.