



ACROS ORGANICS part of Thermo Fisher Scientific

ThermoFisher
S C I E N T I F I C

Product Specifications

Product 41709
Petrolatum, white,pure

General Product Data

Version	0
CAS No	8009-03-8
Molecular weight	0.00
Molecular formula	
Linear formula	
Flash point (°C)	>150

Product Specifications

Appearance (Color)	White
Appearance (Form)	Paste
Infrared spectrum	Conforms
Solidification point	50°C to 56°C
Viscosity	5 to 8 mm ² /s (at 100°C)

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticaaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.thermofisher.com>
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329