

PDS No. 2141x1	PRODUCT DATA SHEET				Page 1 of 1
Revision 01	Container, PS, 31 x 56 mm				 greiner bio-one
	Greiner Item-No. 2141x1				
Valid for Item-No.:	214101	214111			

1.	Description / Specification	
1.1	Description	PS Container, 31 x 56 mm 214101: container without stopper (appropriate multi-disk stopper, Cat.-No. 328312) 214111: container with multi-disk stopper
1.2	Dimensions	<u>Container</u> : height: 56 mm inner diameter: 29 mm outside diameter: 31 mm <u>Stopper</u> : diameter (thread): 29,8 mm (± 0,2 mm) outside diameter: 32,2 mm (± 0,15 mm) Complete (container + stopper): height: 59,2 mm
1.3	Volume	Total volume: 34 ml Working volume: 29 ml
1.4	Material / Resin	<u>Container</u> : PS (Polystyrene), free of heavy metal <u>Stopper</u> : LDPE (Low-density polyethylene), free of heavy metal
1.5	Colour	<u>Container</u> : clear <u>Stopper</u> : white
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	<u>Container</u> : -20°C to +60°C <u>Stopper</u> : -50°C to +80°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	214101	214111
3.1	Pieces / Bag	1300	50
3.2	Pieces / Box	1300	800
3.3	Lot-No.	JJ KW 01 (year, week, machine number)	JJ KW 01 (year, week, machine number)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision -	Date 12 October 2012	Date 17 October 2012	Date 17 October 2012	
Date -	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	