



KLINIPATH Microscope Slides

Description:	Clear imprinted with KLINIPATH LOGO ; Clear with a “●” at each of the two lower corners; pre-cleaned, clear; no marks, prints or contamination, etc. Color: Standard colors Color Frosted Area: 19.5 x 25 (±0.5)mm;
Purpose:	Used for histological pathology examination
Materials:	Standard Glass; thickness 1.0 (-0.05,+0.1)mm
Features:	The smooth edges provide the safety of the operation
	Color Frosted
	Suitable for marking with conventional labeling system, pencil, mark pens and inkjet printers
Temperature:	Normal laboratory temperature and humidity
Uses:	Recommended for routine H&E stains,frozen sections
Shelf life:	60 Months

KLINIPATH Microscope Slides

EU Part No.	Dimension				Frosted or Painted	Paint Color	Corner Angle	Edges	Packaging	QTY/CASE
	A (mm)	B (mm)	C (mm) - 0.05/+0.1	D (mm) (+/-0.5)						
PR-I-001	75.7~76.0	25.0~25.3	1	N/A	Painted	WHITE	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500
PR-I-002	75.7~76.0	25.0~25.3	1	N/A	Painted	BLUE	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500
PR-I-003	75.7~76.0	25.0~25.3	1	N/A	Painted	GREEN	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500
PR-I-004	75.7~76.0	25.0~25.3	1	N/A	Painted	YELLOW	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500
PR-I-005	75.7~76.0	25.0~25.3	1	N/A	Painted	YELLOW	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500
PR-I-006	75.7~76.0	25.0~25.3	1	N/A	Painted	YELLOW	45° (small corner)	Smooth	50pcs / Box, 50 Boxes / Case	2500