





TUBE STYLE SELECTOR

Example Tube Style									
Example Tube/Bottle/Vial Type	Blood Tubes Culture Tubes Test Tubes	Blood Tubes Culture Tubes Test Tubes	Test Tubes Culture Tubes 15mL Tube	Test Tubes Culture Tubes 15mL Tube	Test Tubes Culture Tubes	Test Tubes Culture Tubes	Test Tubes Culture Tubes 50mL Tube	Test Tubes Culture Tubes	
Well Tolerance	10-13 mm	13-16mm	16-18mm	16-20mm	20-24mm	20-25mm	25-30mm	30-40mm	
Rack Codes	HS120089 120759	120761 HS23072 120763 120760	HS23215	120764 120765	120766	120767 120768 120769	120770 HS232461	120771 120772	

Example Tube Style										
Example Tube/Bottle/Vial Type	Bijou Bottle Scintillation Vial		McCarthy Bottle Universal Bottle Scintillation Vial Culture Bottle Serum Vial		Urine Containers Centrifuge Bottles Culture Bottles Storage Bottles			Centrifuge Bottles Culture Bottles Storage Bottles		
Well Tolerance	23-25mm		30-33mm		55-58mm			55-60mm		
Rack Codes	HS120085 HS120086		HS120087 HS120088		HS120091			HS23224		

TUBE ARRAY EXAMPLES

2x3	2x4	2x6	3x5	3x8	4x6	4x10
						
4x12	5x10	6x8	6x12	8x10	9x12	
	