



FP4

Flocculators

Portable digital Jar Tester with brushless motor,
timer and lighting system



Practical and reliable solution for on-site Jar Test

Digital portable flocculators in 4 positions with brushless motor, timer and central lighting system to allow standard conditions to be met during Jar Testing in water treatment plants.

- Identical conditions for all the 4 positions for comparable results
- Digitally selectable speed from 10 to 300 rpm
- Settable timer for unattended operation
- Stirring rods adjustable in height without tools
- Powerful brushless motor and SpeedServo technology
- Central light to simplify sample observation
- Practical handle for on-site testing at water treatment plants
- Strong resistance to chemical and mechanical corrosion



INSTRUMENT - CODE

FP4

100-240 V / 50-60 Hz

F105A0108

TECHNICAL DATA

	FP4
CONSTRUCTION MATERIAL	Epoxy painted metal structure
MOTOR	Brushless
USER INTERFACE	Digital
NUMBER OF SAMPLES	4 positions
LIGHTING SYSTEM	Central light
STAINLESS STEEL STIRRING RODS	Adjustable in height
MONOSELECTOR	Same speed for each rod
ELECTRONIC SPEED CONTROL	From 10 to 300 rpm
SPEED SETTING INTERVAL	1 rpm
TIMER	1÷999 min or 1÷99 hours
POWER	6 W
DIMENSIONS (WxHxD)	250x310x265 mm (9.8x12.2x10.4 in)
WEIGHT	3.8 kg (8.3 lb)



We reserve the right to make technical alterations. We do not assume liability for errors in printing, typing or transmission.

VELP Official Partner



ITALY - HQ

Via Stazione 16 20865
Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA

velpindia@velp.com

USA

155 Keyland Court, Bohemia
NY 11716 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

LATAM

velplatam@velp.com

CHINA

Xinlong Rd Building 28, Lane
1333 Shanghai city - China
Tel. +86 18017557329
velpchina@velp.com