


L06841 2,6-Diisopropylphenol, 97%

Product Number: L06841
CAS number: 2078-54-8
MDL number: MFCD00008885
Molecular formula: $C_{12}H_{18}O$
Molecular weight: 178.28



Product Specification

Appearance (Color): Clear, pale yellow to orange-brown
Form: liquid
Assay (GC): > 96.0%
Refractive Index: 1.5115 - 1.5165 @ 20 °C

Date of Print: June 10, 2021
Version: 8

Product Specifications are subject to amendment and may change over time.