



ACROS ORGANICS part of Thermo Fisher Scientific

ACR^{OS}
ORGANICS

Product Specifications

Product 116220025
N,N-Dimethylformamide, 99+%, extra pure

General Product Data

Version	5
CAS No	68-12-2
Molecular weight	73.09
Molecular formula	C ₃ H ₇ N O
Linear formula	HCON(CH ₃) ₂
Flash point (°C)	58

Product Specifications

Appearance (Form)	Clear liquid
Color scale	=<20 APHA
GC	>=99.0 %
Water	=<0.2 % (K.F.)
Residue after evaporation	=<0.005 %

Product Specifications are subject to amendment and may change over time.

Acros Organics

ENA23, zone 1, nr 1350, Janssen Pharmaceutica laan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: <http://www.acros.com>
1 Reagent Lane, Fair Lawn, NJ 07410, USA
Fax 201-796-1329