

Certificate of Conformity

BOUFFANT CAP GREEN 61 CM – BAG 113-8289

CHARACTERISTICS:

Bouffant caps, also called head covers, are made of 12 grams/m² by weight, Non-woven polypropylene material, commonly used in a clean environment. Bouffant caps are manufactured, vacuum-sealed in a stringent clean environment. Color green.

PHYSICAL PROPERTIES

Material	Non-Woven Polypropylene	
Elastic	latex free	
Size	61 cm (24 inch) + 1.0"	
Stretch Length	61 cm (24 inch) + 1.0"	
Material Thickness	2.0 mils + 0.2 mils	
Colors	Green	
Tensile Strength	MD > 16N/5cm;	CD > 10N/5cm
Ultimate Elongation	>30 %	

PACKING:

100pcs/bag, 10bags/case

QUALITY ASSURANCE CRITERIA:

SPP overalls are guaranteed to meet 4.0 AQL level II sampling plan per MIL-std-150E for the following defect criteria:

- | | |
|--------------------------|--|
| 1. Broken elastic band | 4. Holes/tears in PP material caused by sewing machine |
| 2. Irregular PP material | 5. Physical dimensions non-com |
| 3. Dirty PP material | |

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **November 2015**