

Material No.: 16004-390

Material Description: VWR Microscope Slide, Unimark, 75x25x1mm, Pink, 72/box, 20 boxes/cs

Lot No.: 55160

Compliance and Procedure Result

Product Specifications

- VWR VistaVision Microscope Slide, plain slide
- 75 x 25mm with thickness 1mm
- Country of origin: Germany
- Soda Lime Glass
- pre-cleaned, Pink marking area of approx. 20mm on one side
- Autoclavable
- ready for use

Accuracy

- Tolerance of dimensions: +/- 0.2 mm
- Tolerance of thickness: 0.95-1.05 mm

Compliance

Standard: DIN ISO 8037-1

Signed on behalf of VWR By: Melissa Bigler

Quality Employee

5/18/2026