

Material Number: 20873.180
Material description: Glutaraldehyde 25% in aqueous solution
Grade: GPR RECTAPUR®
CAS Number: 111-30-8
Molecular weight: 100.12 g/mol
Molecular formula: C₅H₈O₂
Date specification set: 3rd Apr 2019

(All batches of this product after this date are tested to this specification. In the unlikely event that the specification is changed a new certificate will be issued)

Characteristics	Specifications
Appearance	clear colorless or yellow liquid
Assay	23 - 27 %
Color scale	max. 50 APHA
Density (d20°C/4°C)	1.055 – 1.065
pH	3.2 – 4.2
Methanol	max. 0,3%

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

VWR International