

Product Specification

Test	Specification
Ammonium (NH ₄)	<= 0.3 ppm
Conductivity (uS/cm)	<= 5
Heavy Metals (as Pb)	passes test
pH at 25°C	5.0 – 7.0
Residue after Evaporation	<= 10 ppm
Substances Reducing KMnO ₄	passes test
Trace Impurities	
Chloride (Cl)	<= 1 ppm
Sulphate (SO ₄)	passes test

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000