

Buffer Solution
Baker Analyzed
pH 7.00 [Phosphate]



Material No: 5608.0500
Effective from: 2025/01/28

Product Specification

Test	Specification
pH at 25°C	6.96 – 7.04
Appearance	passes test

For Laboratory, Research or Manufacturing use.

This buffer is traceable to NIST Standard Reference Buffer. Contains Potassium Dihydrogen Phosphate / Disodium Hydrogen Phosphate. Contains 0.01% Sodium Benzoate.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000