



Material	190063U
Material description	Ammonia 10%
Grade	AnalaR NORMAPUR
CAS Number	1336-21-6
Molecular formula	NH ₄ OH
Molecular mass	35.05

Characteristics	Specifications
Assay	Min. 10.0 %
Appearance	Clear colourless liquid
Cl (Chloride)	Max. 0.5 ppm
PO ₄ (Phosphate)	Max. 0.1 ppm
S (Sulphide)	Max. 0.1 ppm
SO ₄ (Sulphate)	Max. 0.1 ppm
Fe (Iron)	Max. 0.1 ppm
K (Potassium)	Max. 1 ppm
Mg (Magnesium)	Max. 0.5 ppm
Mn (Manganese)	Max. 0.1 ppm
Na (Sodium)	Max. 1 ppm
Ni (Nickel)	Max. 0.05 ppm
Pb (Lead)	Max. 0.05 ppm
Zn (Zinc)	Max. 0.2 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.