



Product Specification Sheet

VWR MICROCENTRIFUGE TUBES WITH HINGED LID

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity	Type	Colour	Packed	Pk
525-0660	5.0 ml	non sterile	Clear	200/pack	200
525-0661	5.0 ml	non sterile	Blue	200/pack	200
525-0662	5.0 ml	non sterile	Yellow	200/pack	200
525-0663	5.0 ml	non sterile	Green	200/pack	200
525-0664	5.0 ml	non sterile	Rose	200/pack	200
525-0665	5.0 ml	non sterile	Black	200/pack	200
525-0666	5.0 ml	sterile	Natural	25/pack	200

Material: Polypropylene - copolymer

Specifications:

Clean Room: Class 100,000

Autoclavable: 120°C (15 minutes)

Centrifugation: Only suitable for centrifuges with swing-out rotor and adapters for 15 ml conical tubes

Max. RCF : 6000 x g

Freezing: -90°C

Sterile (Optional): Gamma-Irradiation

Other Characteristics:

- Good chemical resistance
- Frosted marking area on the tube cap is provided for sample tagging and identification
- RNase & DNase Free / Pyrogen free for the sterile version

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **October 2014**